

Comfortland Ankle Gauntlet

BEFORE USING THIS DEVICE, PLEASE READ THE FOLLOWING INSTRUCTIONS COMPLETELY AND CAREFULLY. CORRECT APPLICATION IS VITAL TO THE PROPER FUNCTIONING OF THIS DEVICE.

INDICATIONS


The Comfortland Ankle Gauntlet provides compression and insulation for diabetics suffering from poor blood circulation and pain in the feet. Effective in raising foot temperatures by 2-3 degrees in an effort to relieve pain.

CAUTION

The application instructions listed above are intended for the use of qualified professionals only. It is not recommended for patients to attempt to apply or adjust this device unless otherwise instructed to do so by an orthopedic professional. Any questions regarding further adjustment or frequency/ duration of wear should be directed to your orthopedic professional, as they know the individual condition for which this device was prescribed.

APPLICATION

1. Undo the closure strap and open the back of the ankle gauntlet. Place your foot in the gauntlet the same as you would put on a sock (Fig. A).
2. Wrap the closure strap around the back of the leg and attach the strap to itself using the Velcro tab on the end of the strap (Fig. B).
3. Tighten and adjust as needed. If properly applied, the Ankle Gauntlet will fit snug and should resemble Fig. C



CLEANING INSTRUCTIONS

Hand wash in cold water with mild detergent. Allow all parts to air dry before re-assembly. Do not dry clean, iron, or bleach any part of this brace, as it could cause damage and void the manufacturer's warranty.

NOTICE: While every effort has been made to obtain maximum strength, function, and durability, there is no guarantee that injury will be prevented through the use of this device.